

# siPOOLS: Complex and Defined siRNA Pools for Specific and Potent Gene Knock-down


Catherine Goh<sup>1</sup>, Michaela Beitzinger<sup>1</sup>, Andrew Walsh<sup>1</sup>, Gunter Meister<sup>2</sup>, Stefan Hannus<sup>1</sup>, Michael Hannus<sup>1</sup>

<sup>1</sup>siTOOLS Biotech GmbH, Lochhamerstrasse 29A, 82152 Martinsried/Planegg, Germany; <sup>2</sup>Biochemistry Center Regensburg (BZR), Laboratory for RNA Biology, University of Regensburg, 93053 Regensburg, Germany.

RNAi is an excellent method for gene function analysis: it is reversible, dose-dependent, versatile, fast and cost-efficient. These drug-like properties have made it the tool of choice for functional genomic screens. However, off-target effects dominate published RNAi screens and limit their yield to a small number of strongly responding, frequently well-established gene factors. siPOOLS are complex pools of 30 carefully selected siRNAs. The very low concentration of each individual siRNA dilutes off-target effects below detection. Target gene silencing becomes highly efficient and reliable due to the cooperative activity of the pooled siRNAs. siPOOLS therefore allow a far deeper interpretation of RNAi screening data, increasing the yield of novel targets and aiding the understanding of complex biology.

## siRNAs: Noise from off-target effects

### Correlation Analysis of Typical RNAi Screening Data




Single dot = **two replicate** values for **one siRNA**

- High correlation
- Great assay reproducibility

Single dot = values of **two siRNAs** targeting **one gene**

- No correlation
- Non-interpretable results

### RNAi Screening is Highly Inefficient




- Two siRNA libraries were screened in an identical assay by the same lab.
- Overlap between hits at noise level (16 expected by chance).
- Majority of hits expected to be false.

- Current RNAi reagents are plagued by off-target effects.
- RNAi screening with current reagents is highly inefficient and unreliable.

## siPOOLS: Complex and defined siRNA pools


### The siPOOL Concept



- An **siRNA** produces many off-target effects when administered alone ⇒ low specificity, unreliable results.
- A **siPOOL** is 30 selected siRNAs, each siRNA at picomolar concentrations ⇒ Minimal off-target effects, co-operative target knock-down, enhanced reliability.

### Why 30 siRNAs?

#### Off-target spiking experiment



- An **siRNA** with known off-target gene (MAD2) was spiked into pools of increasing complexity.
- Luciferase assay detects activity of off-target gene.
- Efficient off-target dilution requires pools of > 15 siRNAs.


- Complex pools of more than 15 siRNAs are required to prevent off-target effects.

... because 4 are not enough

## siPOOLS: Specific and potent gene knock-down

### Specificity


#### Whole genome expression analysis



- An **siRNA** (left) induces massive off-target gene deregulation.
- The same siRNA incorporated into a **siPOOL** (right) shows dramatically reduced off-targets.

### Potency


#### siPOOL-mediated gene knockdown: 1 nM is sufficient



- Most **siPOOLS** achieve 75-98% target gene knockdown at less than 1 nM.

### Reproducible Phenotypes

#### Phenotypic correlation study



- 36 genes were silenced with **3 siRNAs/gene** (left) or **2 siPOOLS/gene** (right).
- Cell viability was measured after 72h.
- Phenotypic correlation was greatly improved with siPOOLS.

### References

- Marine S, Bahl A, Ferrer M, Buehler E (2012). Common seed analysis to identify off-target effects in siRNA screens. *Journal of Biomolecular Screening*, 17(3), 370-8.
- Hasson SA, Kane LA, Yamano K, Huang CH, Sliter DA, Buehler E, Wang C, Heman-Ackah SM, Hessa T, Guha R, Martin SE, Youle RJ (2013). High content genome-wide RNAi screens identify regulators of parkin upstream of mitophagy. *Nature*, 504(7479), 291-5.
- Hannus M, Beitzinger M, Engelmann JC, Weickert MT, Spang R, Hannus S, Meister G (2014). siPools: highly complex but accurately defined siRNA pools eliminate off-target effects. *Nucleic Acids Research*, 42(12), 8049-61.

- Specific and potent gene knock-down with siPOOLS ensures reproducible and reliable phenotypes.
- RNAi screening with siPOOLS produces high-quality hits.