

# siPOOLS: highly complex but accurately defined siRNA pools eliminate off-target effects

**Abstract:** Short interfering RNAs (siRNA) are widely used as tool for gene inactivation in basic research and therapeutic applications. One of the major shortcomings of siRNA experiments are sequence-specific off-target effects. Such effects are largely unpredictable because siRNAs can affect partially complementary sequences and function like microRNAs (miRNAs), which inhibit gene expression on mRNA stability or translational levels. Here we demonstrate that novel, complex siRNA pools - referred to as siPOOLS -


containing up to 60 selected siRNAs eliminate off-target effects while maintaining highly efficient and reliable target gene silencing. This is achieved by the low concentration of each individual siRNA diluting sequence-specific off-target effects below detection limits. In fact, whole transcriptome analysis reveal that single siRNA transfections can severely affect global gene expression. However, when complex siRNA pools are transfected, almost no transcriptome alterations are observed.

Josef Unger<sup>a</sup>, Michael Hannus<sup>a,b</sup>, Michaela Beitzinger<sup>a,b</sup>, Julia C. Engelmann<sup>c</sup>, Marie-Theresa Weickert<sup>b</sup>, Rainer Spang<sup>d</sup>, Stefan Hannus<sup>c</sup> and Gunter Meister<sup>b</sup>

a: siTOOLS Biotech GmbH, Lochhamerstrasse 29A, 82152 Planegg;  
b: Biochemistry Center Regensburg (BZR), Laboratory for RNA Biology, University of Regensburg, Universitätsstraße 31, 93053 Regensburg;


c: Intana Bioscience GmbH, Lochhamerstraße 29A, 82152 Planegg;  
d: Department of Statistical Bioinformatics, Institute for Functional Genomics; University of Regensburg, Josef-Engert-Straße 9, 93053 Regensburg

## Correlation analysis of siRNA screening results




➔ Single siRNA Phenotypes are reproducible but are dominated by off-target effects

Reference: Common Seed Analysis to Identify Off-Target Effects in siRNA Screens  
*J Biomol Screen* 2011, Shane Marine, Amit Bahl, Marc Ferrer and Eugene Buehler



## siRNAs trigger off-target effects acting like microRNAs



## >15 siRNAs are required for efficient off-target dilution


1: siRNAs with known strong off-target (King et al. 2012) for pool spiking




2: Luciferase off-target effect reporter




3: Results



## siPools: Concept and Technology



## siPOOLS do not affect global gene expression



Affymetrix array: 24h post transfection, 10nM, HeLa