

Jun 14, 2016





Welcome to our first siTOOLS newsletter: Stay updated with the latest company and product-related news.

siPOOLS/siPOOL libraries for high specificity RNAi validation/screening

Beyond individual siPOOLS, the first siPOOL library is now available targeting 505 classical human protein kinases. Choose the amount you need and define the layout that works best for your screening assay. Customized siPOOL libraries are also available.


siPOOLS give reliable phenotypes – Test screen with 36 genes from the kinase library:



If reliable and specific, two reagents that target the same gene should produce similar phenotypes. The siPOOLS clearly outperformed individual siRNAs in achieving similar phenotypes.

Request a quote or find out more!

New Product - raPOOLS for RNA affinity purification



Discover interacting proteins, RNA or DNA - affordable and easy RNA pulldown.

An RNA antisense pool (raPOOL) combines 30 biotinylated probes, synthesized and designed by proprietary means, enabling specific and robust enrichment of target RNA at non-prohibitive cost. Save time designing probes and achieve unparalleled value with prices as low as 8 Euro/probe.

Request a quote or find out more!

Pricing and Ordering

We have recently simplified our pricing and ordering formats, making siPOOLS even more affordable compared to conventional RNAi reagents. Now, you can assemble your own siPOOL kit and have more siPOOLS at a discounted price.

To take advantage, visit our webshop!

Publications and Testimonials

siPOOLS are used by a growing number of scientists in academic and commercial biomedical research. [Here](#) you will find recent publications using siPOOLS. Also check out our latest poster on RNAi and CRISPR. Have we missed your publication involving siPOOLS? Please let us know and we will be happy to highlight your paper on our website!

We are also collecting testimonials that will be featured on our website. Help us spread the word so that others like you can enjoy the benefit and utility of working with siPOOLS.

Send your testimonial now!



©siTOOLS Biotech GmbH. All rights reserved.

www.sitoolsbiotech.com
blog.sitoolsbiotech.com

info@sitoolsbiotech.com

+49 (0) 89 4431 2584